

ET Series - Ratio Measurement Filter Sets




2023 V1


For customized projects please Contact us:

[**info@simtrum.com**](mailto:info@simtrum.com)

79001




79002




FILTERS	TYPE	%T-OD	AOI	DATA	FILTERS	TYPE	%T-OD	AOI	DATA
● ET340x	EX	0.00 T at 300nm	0°	ASCII	● ET340x	EX	0.00 T at 300nm	0°	ASCII
● ET380x	EX	0.00 T at 300nm	0°	ASCII	● ET380x	EX	0.00 T at 300nm	0°	ASCII
● T400lp	BS	0.00 T at 300nm	45°	ASCII	● ET436/20x	EX	0.00 T at 300nm	0°	ASCII
● ET510/80m	EM	0.00 T at 300nm	0°	ASCII	● ET495/10x	EX	0.00 T at 300nm	0°	ASCII
					● T510lpxr	BS	0.00 T at 300nm	45°	ASCII
					● ET560/80m	EM	0.00 T at 300nm	0°	ASCII

79003




79004




FILTERS	TYPE	%T-OD	AOI	DATA	FILTERS	TYPE	%T-OD	AOI	DATA
● ET340x	EX	0.00 T at 300nm	0°	ASCII	● ET340x	EX	0.00 T at 300nm	0°	ASCII
● ET380x	EX	0.00 T at 300nm	0°	ASCII	● ET380x	EX	0.00 T at 300nm	0°	ASCII
● ET555/25x	EX	0.00 T at 300nm	0°	ASCII	● ET470/40x	EX	0.00 T at 300nm	0°	ASCII
● 79003bs	BS	0.00 T at 300nm	45°	ASCII	● T495lpxr	BS	0.00 T at 300nm	45°	ASCII
● 79003m	EM	0.00 T at 300nm	0°	ASCII	● ET535/70m	EM	0.00 T at 300nm	0°	ASCII

79005




FILTERS	TYPE	%T-OD	AOI	DATA
● ET436/20x	EX	0.00 T at 300nm	0°	ASCII
● ET495/10x	EX	0.00 T at 300nm	0°	ASCII
● T510lpxru	BS	0.00 T at 300nm	45°	ASCII
● ET540/40m	EM	0.00 T at 300nm	0°	ASCII

79006




FILTERS	TYPE	%T-OD	AOI	DATA
● ET365/10x	EX	0.00 T at 300nm	0°	ASCII
● T387lp	BS	0.00 T at 300nm	45°	ASCII
● ET405/30m	EM	0.00 T at 300nm	0°	ASCII
● ET485/25m	EM	0.00 T at 300nm	0°	ASCII

79009



FILTERS	TYPE	%T-OD	AOI	DATA
● ET500/20x	EX	0.00 T at 450nm	0°	ASCII
● T570lpxr	BS	0.00 T at 450nm	45°	ASCII
● ET585/20m	EM	0.00 T at 450nm	0°	ASCII
● ET640/20m	EM	0.00 T at 450nm	0°	ASCII

79010



FILTERS	TYPE	%T-OD	AOI	DATA
● ET550/20x	EX	0.00 T at 500nm	0°	ASCII
● T570lpxr	BS	0.00 T at 500nm	45°	ASCII
● ET585/20m	EM	0.00 T at 500nm	0°	ASCII
● ET640/20m	EM	0.00 T at 500nm	0°	ASCII