

Laser Bandpass Filter Set




2023 V1

For customized projects please Contact us:

info@simtrum.com


ET Series - Laser Single Band Bandpass Filter Sets

49901




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/20x	EX	0.00 T at 300nm	0°	ASCII
● ZT405rdc	BS	0.00 T at 300nm	45°	ASCII
● ET460/50m	EM	0.00 T at 300nm	0°	ASCII

49902




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET436/20x	EX	0.00 T at 350nm	0°	ASCII
● ZT442rdc	BS	0.00 T at 350nm	45°	ASCII
● ET480/40m	EM	0.00 T at 350nm	0°	ASCII

49903




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET436/20x	EX	0.00 T at 350nm	0°	ASCII
● ZT442rdc	BS	0.00 T at 350nm	45°	ASCII
● ET500/40m	EM	0.00 T at 350nm	0°	ASCII

49904




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/10x	EX	0.00 T at 400nm	0°	ASCII
● ZT488rdc	BS	0.00 T at 400nm	45°	ASCII
● ET525/50m	EM	0.00 T at 400nm	0°	ASCII

49905



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET514/10x	EX	0.00 T at 450nm	0°	ASCII
● ZT514rdc	BS	0.00 T at 450nm	45°	ASCII
● ET545/40m	EM	0.00 T at 450nm	0°	ASCII


49906



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET514/10x	EX	0.00 T at 450nm	0°	ASCII
● ZT514rdc	BS	0.00 T at 450nm	45°	ASCII
● ET575/50m	EM	0.00 T at 450nm	0°	ASCII


ET Series - Laser Single Band Bandpass Filter Sets

49907




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET532/10x	EX	0.00 T at 450nm	0°	ASCII
● ZT532rdc	BS	0.00 T at 450nm	45°	ASCII
● ET575/50m	EM	0.00 T at 450nm	0°	ASCII

49908




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET532/10x	EX	0.00 T at 450nm	0°	ASCII
● ZT532rdc	BS	0.00 T at 450nm	45°	ASCII
● ET590/50m	EM	0.00 T at 450nm	0°	ASCII

49909




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET561/10x	EX	0.00 T at 500nm	0°	ASCII
● ZT561rdc	BS	0.00 T at 500nm	45°	ASCII
● ET600/50m	EM	0.00 T at 500nm	0°	ASCII

49910




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET561/10x	EX	0.00 T at 500nm	0°	ASCII
● ZT561rdc	BS	0.00 T at 500nm	45°	ASCII
● ET630/75m	EM	0.00 T at 500nm	0°	ASCII

49911



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET594/10x	EX	0.00 T at 500nm	0°	ASCII
● ZT594rdc	BS	0.00 T at 500nm	45°	ASCII
● ET645/75m	EM	0.00 T at 500nm	0°	ASCII


49914



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET642/20x	EX	0.00 T at 600nm	0°	ASCII
● ZT647rdc	BS	0.00 T at 600nm	45°	ASCII
● ET705/72m	EM	0.00 T at 600nm	0°	ASCII


ET Series - Laser single-band long-pass filter sets

49912




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET594/10x	EX	0.00 T at 550nm	0°	ASCII
● ZT594rdc	BS	0.00 T at 550nm	45°	ASCII
● ET610lp	EM	0.00 T at 550nm	0°	ASCII

49913



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET635/20x	EX	0.00 T at 550nm	0°	ASCII
● ZT640rdc	BS	0.00 T at 550nm	45°	ASCII
● ET655lp	EM	0.00 T at 550nm	0°	ASCII


49915



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET365/20x	EX	0.00 T at 300nm	0°	ASCII
● ZT365bcm	BS	0.00 T at 300nm	45°	ASCII
● ET425lp	EM	0.00 T at 300nm	0°	ASCII


ET Series - Laser Multi-Band Filter Sets

59901




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/488x	EX	0.00 T at 350nm	0°	ASCII
● ZT405/488rpc	BS	0.00 T at 350nm	45°	ASCII
● ZET405/488m	EM	0.00 T at 350nm	0°	ASCII

59902




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET442/514x	EX	0.00 T at 400nm	0°	ASCII
● ZT442/514rpc	BS	0.00 T at 400nm	45°	ASCII
● ZET442/514m	EM	0.00 T at 400nm	0°	ASCII

59903




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET442/532x	EX	0.00 T at 400nm	0°	ASCII
● ZT442/532rpc	BS	0.00 T at 400nm	45°	ASCII
● ZET442/532m	EM	0.00 T at 400nm	0°	ASCII

59904




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/561x	EX	0.00 T at 450nm	0°	ASCII
● ZT488/561rpc	BS	0.00 T at 450nm	45°	ASCII
● ZET488/561m	EM	0.00 T at 450nm	0°	ASCII

59905



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/594x	EX	0.00 T at 450nm	0°	ASCII
● ZT488/594rpc	BS	0.00 T at 450nm	45°	ASCII
● ZET488/594m	EM	0.00 T at 450nm	0°	ASCII


59906



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/640x	EX	0.00 T at 450nm	0°	ASCII
● ZT488/640rpc	BS	0.00 T at 450nm	45°	ASCII
● ZET488/640m	EM	0.00 T at 450nm	0°	ASCII


ET Series - Laser Multi-Band Filter Sets

59901




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/488x	EX	0.00 T at 350nm	0°	ASCII
● ZT405/488rpc	BS	0.00 T at 350nm	45°	ASCII
● ZET405/488m	EM	0.00 T at 350nm	0°	ASCII

59902




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET442/514x	EX	0.00 T at 400nm	0°	ASCII
● ZT442/514rpc	BS	0.00 T at 400nm	45°	ASCII
● ZET442/514m	EM	0.00 T at 400nm	0°	ASCII

59903




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET442/532x	EX	0.00 T at 400nm	0°	ASCII
● ZT442/532rpc	BS	0.00 T at 400nm	45°	ASCII
● ZET442/532m	EM	0.00 T at 400nm	0°	ASCII

59904




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/561x	EX	0.00 T at 450nm	0°	ASCII
● ZT488/561rpc	BS	0.00 T at 450nm	45°	ASCII
● ZET488/561m	EM	0.00 T at 450nm	0°	ASCII

59905



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/594x	EX	0.00 T at 450nm	0°	ASCII
● ZT488/594rpc	BS	0.00 T at 450nm	45°	ASCII
● ZET488/594m	EM	0.00 T at 450nm	0°	ASCII


59906



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET488/640x	EX	0.00 T at 450nm	0°	ASCII
● ZT488/640rpc	BS	0.00 T at 450nm	45°	ASCII
● ZET488/640m	EM	0.00 T at 450nm	0°	ASCII


ET Series - Laser Multi-Band Filter Sets

59907




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET532/640x	EX	0.00 T at 500nm	0°	ASCII
● ZT532/640rpc	BS	0.00 T at 500nm	45°	ASCII
● ZET532/640m	EM	0.00 T at 500nm	0°	ASCII

69901




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/488/561x	EX	0.00 T at 350nm	0°	ASCII
● ZT405/488/561rpc	BS	0.00 T at 350nm	45°	ASCII
● ZET405/488/561m	EM	0.00 T at 350nm	0°	ASCII

69902




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/488/594x	EX	0.00 T at 350nm	0°	ASCII
● ZT405/488/594rpc	BS	0.00 T at 350nm	45°	ASCII
● ZET405/488/594m	EM	0.00 T at 350nm	0°	ASCII

69904




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET442/514/561x	EX	0.00 T at 400nm	0°	ASCII
● ZT442/514/561rpc	BS	0.00 T at 400nm	45°	ASCII
● ZET442/514/561m	EM	0.00 T at 400nm	0°	ASCII

69905



FILTERS	TYPE	%T-OD	AOI	DATA
● ZET445/514/594x	EX	0.00 T at 351nm	0°	ASCII
● ZT445/514/594rpc	BS	0.00 T at 351nm	45°	ASCII
● ZET445/514/594m	EM	0.00 T at 351nm	0°	ASCII


89901v2




FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/488/561/640xv2	EX	0.00 T at 350nm	0°	ASCII
● ZT405/488/561/640rpcv2	BS	0.00 T at 350nm	45°	ASCII
● ZET405/488/561/640mv2	EM	0.00 T at 350nm	0°	ASCII

ET Series - Laser Multi-Band Filter Sets

89902



89904



FILTERS	TYPE	%T-OD	AOI	DATA	FILTERS	TYPE	%T-OD	AOI	DATA
● ZET405/488/561/647x	EX	0.00 T at 350nm	0°	ASCII	● ZET405/445/514/561/640x	EX	0.00 T at 350nm	0°	ASCII
● ZT405/488/561/647rpc	BS	0.00 T at 350nm	45°	ASCII	● 89903bs	BS	0.00 T at 350nm	45°	ASCII
● ZET405/488/561/647m	EM	0.00 T at 350nm	0°	ASCII	● ZET405/445/514/561/640m	EM	0.00 T at 350nm	0°	ASCII